

10 mm

ET file system

TransformX
root canal shaping system

Technique Workflow

10 mm

30 mm

44mm

STEP 1

STEP 2

STEP 3

30 mm



77 mm

57 mm

TransformX™
Acrobat
13/.03

GLIDE PATH

300 RPM
2.0~2.5 Ncm

TransformX™
ET 20/04

SHAPING

350 RPM
3.0 Ncm

TransformX™
ET 25/.04

FINISHING

350 RPM
3.0 Ncm
(Option: ET runs
up to 500 RPM)

Irrigate → Recapitulate → Re-Irrigate

10 mm

10

30 mm

30 mm